

## BENCHMIXER™ XL Multi-Tube Vortexer!




† BenchMixer™ XL은 tube, flask, 시약병, 실린더에 담긴 시료를 핸드 프리로 믹싱 작업에 적합.

**QuEChERS 잔류농약등 유해물질 분석 전처리에 적합하며,**

매우 다양한 종류의 용기, rack, 기타 응용 기구 (수평 또는 수직)에 바로 적용가능. 펄싱 프로그램 기능은 시간 간격을 두고 특정 시간에만 작동 하면서 교반 동작의 강도 높일수 있다. 버튼만 누르면 최고 100시간까지 작업자 감독 없이 동작.



50 X 1.5/2.0ml



15 X 50ml



50 X 5.7 or 10ml



3 x microplates/96 Well Tube Racks



6 X 50, 100 or 250ml Volumetric Flasks



9 X 50ml Horizontal QuEChERS Method

### Technical Data

Speed range :	500 to 2500rpm
Max. capacity:	6 x 250ml
Orbit:	3.6mm
Display:	Timer, rpm digital display
Timer:	1min. to 99 hour 59min. / continuous
Platform size:	30 x 17.8 cm
Dimension (w x d x h)cm:	42 x 25 x 48
Weight 및 Electrical:	14kg / 115V or 230V, 50/60 Hz, 60W

### Ordering Information

BV1010	BenchMixer™ XL Multi-Tube Vortexer 및 50x12mm tube rack 기본포함
BV1010-00	Replacement tray pad set
BV1010-05	Tube rack, 72 x 0.5ml, (Ø: 8mm)
BV1010-1520	Tube rack, 50 x 1.5ml/2.0, (Ø: 10mm)
BV1010-12	Tube rack, 50 x 12mm (standard included)
BV1010-13	Tube rack, 50 x 13mm
BV1010-150	Tube rack, 50 x 15ml, (Ø:16.58mm)
BV1010-500	Tube rack, 15 x 50ml, (Ø: 29mm)
BV1010-25	Tube rack, 15 x 25mm
BV1010-50H	Tube rack, 9 x 50ml 수평형,
BV1010-V50	Rack, 25 x 50ml, volumetric flask
BV1010-V125	Rack, 12x 125ml, volumetric flask
BV1010-V250	Rack, 8 x 250ml, volumetric flask
BV1010-MP	Rack, Max. 3 microplate 또는 96well micro tube rack
BV1010-MP02	Carrier for above, 96 x 0.2ml tubes or Strips (3pk.)